COSMOS PLUS GENGIVA

PORTUGUÊS / ENGLISH / ESPAÑOL

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1.2. Product trade name

- Cosmos Plus Gingiva 500g.

1.3. Description of the physical principle and fundamentals of the product's technology, applied to its operation and action.

Cosmos Plus Gingiva Resin is a photopolymerizable resin, available in pink (medium gingiva), for the manufacture of gingival models, which should be made by professionals qualified in 3D printing using DLP and LCD technologies.

1.4 Technical specifications and characteristics

1.4.1 Composition

Components
Methacrylate oligomers
Methacrylate monomers
Photoinitiator
Pigments

2. STORAGE

Store in a dry place at a temperature between 35.6°F and 86°F. The Cosmos Plus Gingiva product must not be exposed to heat and sunlight, or other sources of direct light.

3. INSTRUCTIONS FOR USE

3.1. Preparing the 3D printing file

- Prepare the STL file in the software related to the 3D printer. Afterwards, it is suggested to initially print a calibration piece available on the Yller website, using the indicated initial parameters.

3.2. Mixing the resin

Before adding the resin to the printer reservoir, shake the resin bottle until it is completely homogenized.

After immersing the resin in a water bath, leave it in the temperature range of 73.4°F to 77°F to facilitate 3D printing.

Add the amount of resin to the reservoir according to its level.

For the resin already in the reservoir, use a silicone spatula and mix it gently to homogenize.

3.3. Post-processing

Once the 3D printing process is complete, remove the printed devices

from the platform. Then carry out the steps described below.

Cleaning: Place the device inside a container and add 92.8% INPM alcohol or isopropyl alcohol until it covers the device and leave for 5 minutes, then remove and place it in a new solution of 92.8% INPM alcohol or isopropyl alcohol for 5 minutes. Afterwards, dry the device with compressed air or let it dry naturally.

- **Pionext Post-Curing:** After drying, post-cure the device in the UV light chamber, using all LEDs on for 5 minutes.

- **Pcure Post-Curing:** After drying, post-cure the device in the UV light chamber, using "Upper" + "Lower" 100% mode for 10 minutes.

- **Anycubic Post-Curing:** After drying, post-cure the device in the UV light chamber, using all LEDs on for 60 minutes. 30 minutes for each side.

- **72W or 36W UV Chamber Post-Curing:** After drying, post-cure the device in the UV light chamber, using all LEDs on for 60 minutes. 30 minutes for each side.

- **Otoflash Post-Curing:** After drying, post-cure the device in the Otoflash, using 2000 cycles on each side, totalling 4000 cycles.

- **Note:** Use the equipment indicated for post-curing or equipment with equivalent specifications.

Finish: Polish to remove support marks.

4. PRECAUTIONS, RESTRICTIONS AND WARNINGS

4.1. Precautions

- When handling this product, it is recommended to use Personal Protective Equipment, such as protective gloves, protective clothing and eye protection.

- Make sure you carry out all the procedures mentioned in the instructions. Inefficient completion of any post-impression step may leave traces of uncured liquid material that

may alter the planning conditions of the device.

- Improper use of the printing equipment may compromise the physical and mechanical properties of the final device produced from this material. Follow the printer's instructions for use provided by the manufacturer and check the need for maintenance and calibration of the equipment.

- Keep the container closed when not in use.

- Keep this product out of the reach of children and animals.

- Do not use the product if its packaging is damaged.

- Do not use the product if it has expired.

- It is the responsibility of the responsible professional to use the products and equipment in accordance with their Instructions for Use.

- Do not use the product without knowledge of the appropriate techniques. Always carry out the correct procedures in the right conditions and with the right equipment.

- Follow the Instructions for Use provided by the manufacturers of the software used for device planning.

- In the event of an accident involving this product, inform the manufacturer using the communication channels available on the website.

- Wash hands and other exposed areas thoroughly after handling.

- Avoid release into the environment.

4.2. Restrictions and Warnings

The use of this product is contraindicated for people with signs of allergy or irritation to the chemical components of the material and/or methacrylates.

5. INDICATION

Product indicated for making gingival models.

6. VALIDITY

Indicated on the packaging.

7. PRODUCT REGISTRATION

Product exempt from ANVISA registration.

8. DISPOSING OF MATERIALS

To dispose of materials, solid waste or remnants of the product, it is

recommended to consult and comply with current legislation. Dispose of liquid contents and container at a special or hazardous waste collection point in accordance with local, regional, national and/or international regulations, at an authorized centre.

If there is no legislation in force, pack them in a sharps collector and dispose of them in hospital waste.

Yller Biomateriais S/A is not responsible for damages, of any nature, caused by incorrect use or use not provided for in the Instructions for Use.

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