COSMOS PLUS CROWN & BRIDGE

PORTUGUÊS / ENGLISH / ESPAÑOL

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- Cosmos Plus Crown & Bridge B1 (50g, 250g, 500g)

- Cosmos Plus Crown & Bridge Bleach (50g, 250g, 500g)

- Cosmos Plus Composite Cleaning Solution.

1.3. Description of the Physical Principle and Fundamentals of the Product's Technology Applied to its Operation and Action

Cosmos Plus Crown & Bridge Micro nano-hybrid is a photopolymerizable resin for the manufacture of permanent bridges, crowns, veneers, inlays and onlays, which must be made by professionals qualified in 3D printing using DLP and LCD technologies.

1.4. Technical specifications and characteristics

1.4.1. Composition

Components
Methacrylate oligomers
Photoinitiator
Inorganic Filler
Additives
Pigments

2. STORAGE

Store in a dry place at a temperature between 35.6°F and 86°F. The product must not be exposed to sunlight.

3.INSTRUCTIONS FOR USE

3.1. Preparing the 3D printing file

- Prepare the STL file in the software related to the 3D printer. Afterwards, it is suggested to initially print a calibration piece available on the Yllér website, using the indicated initial parameters.

- The recommended minimum thickness is 1.4 mm, except for ve-

neers.

- In crowns, do not exceed 3 element units.

3.2. Mixing the resin

- Before adding the resin to the printer's reservoir, place the bottle in a roller agitator for 30 minutes so that the inorganic filler particles are thoroughly distributed.
- After bringing the resin to room temperature in a water bath, keep it within the temperature range of 82.4°F to 89.6°F to facilitate 3D printing.
- Add the appropriate amount of resin to the reservoir, depending on its current level. After printing is complete, return the resin to the bottle.

Note: For the 50g version, add the total amount directly to the reservoir to make the device.

3.3. Post-processing

- Once the 3D printing process is complete, remove the printed devices from the platform and carefully remove the supports. Then use one of the methods described below to remove excess resin from the device.

- **Option 1:** Place the device in a beaker and add Cosmos Plus Composite Cleaning Solution until it covers the device and leave it for 1 minute, then brush it with a toothbrush for 1 minute. This process should be repeated 3 times, always changing the Cosmos Plus Composite Cleaning Solution. Afterwards, dry the device with compressed air.

- **Option 2:** Place the device in a beaker and add Cosmos Plus Composite Cleaning Solution until it covers the device and leave it for 3 minutes in the ultrasound, then remove it from the ultrasound and place it in a new Cosmos Plus Composite Cleaning Solution for 3 minutes at rest. Afterwards, dry the device with compressed air.

- **Pcure Post-Curing:** Immediately after drying with compressed air,

post-cure the device using the "Upper + Lower 100%" mode for 20 minutes.

- Anycubic Post-Curing:

Immediately after drying with compressed air, post-cure the device for 60 minutes.

- **72 W or 36 W Chamber Post-Curing:** Immediately after drying with compressed air, post-cure the device for 60 minutes.

- Otofash Post-Curing:

Immediately after drying with compressed air, post-cure the device using 2000 on each side, totalling 4000 cycles.

Note: Use the equipment indicated for post-curing or equipment with equivalent specifications.

- Finish: Polish to remove the support marks and then carry out the make-up process, according to the manufacturer's instructions.
- Finally, immerse the device in 70° alcohol for 10 minutes.

3.4. Cementation

- Carry out the protocol for preparing the prosthetic pieces.
- Apply Sylano (Yllor) according to the Instructions for Use.
- Prepare the teeth structure to receive the prosthetic piece and carry out the adhesive protocol. For best performance, we recommend using Astra Universal (Yllor). Follow the Instructions for Use.
- Proceed to cement the piece using the Ycem 4All Dual Resin Cement (Yllor) according to the Instructions for Use.

4. PRECAUTIONS, RESTRICTIONS AND WARNINGS

4.1. Precautions

- When handling this product, it is recommended to use Personal Protective Equipment, such as protective gloves, protective clothing, and eye protection.
- Make sure you carry out all the procedures mentioned in the instructions. Inefficient completion of any post-impression step may

leave traces of uncured liquid material that may alter the planning conditions of the device.

- Improper use of the printing equipment may compromise the physical and mechanical properties of the final device produced from this material. Follow the printer's instructions for use provided by the manufacturer and check the need for maintenance and calibration of the equipment.

- Keep the container closed when not in use.

- Keep this product out of the reach of children and animals.

- Do not use the product if its packaging is damaged.

- Do not use the product if it has expired.

- It is the responsibility of the professional responsible to use the products and equipment in accordance with their Instructions for Use.

- Do not use the product without knowledge of the appropriate techniques. Always carry out the correct procedures in the right conditions and with the right equipment.

- Follow the Instructions for Use provided by the manufacturers of the software used for device planning.

- In the event of an accident involving this product, inform the manufacturer using the communication channels available on the website.

- Wash hands and other exposed areas thoroughly after handling.

- Avoid release into the environment.

4.2. Restrictions and warnings

The use of this product is contraindicated for people with signs of allergy or irritation to the chemical components of the material and/or methacrylates.

5. INDICATION

Product indicated for making permanent bridges, crowns, veneers, inlays and onlays.

6. VALIDITY

Indicated on the packaging.

7. PRODUCT REGISTRATION

Reg. Anvisa No: 81009359003.

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- Cosmos Plus Crown & Bridge A1 (50g, 250g, 500g)
- Cosmos Plus Crown & Bridge A2 (50g, 250g, 500g)
- Cosmos Plus Crown & Bridge A3 (50g, 250g, 500g)
- Cosmos Plus Crown & Bridge B1 (50g, 250g, 500g)
- Cosmos Plus Crown & Bridge Bleach (50g, 250g, 500g)
- Cosmos Plus Composite Cleaning Solution.

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