COSMOS PLUS VENEER

PORTUGUÊS / ENGLISH / ESPAÑOL

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1.2 Product trade name

- Cosmos Plus Veneer A1 50g

- Cosmos Plus Veneer A1 250g
- Cosmos Plus Veneer A1 500g
- Cosmos Plus Veneer B1 50g
- Cosmos Plus Veneer B1 250g
- Cosmos Plus Veneer B1 500g
- Cosmos Plus Veneer BL2 50g
- Cosmos Plus Veneer BL2 250g
- Cosmos Plus Veneer BL2 500g

1.3. Description of the physical principle and fundamentals of the product's technology, applied to its operation and action.

Cosmos Plus Veneer is a light-curing, nanohybrid resin for making veneers. They must be manufactured by professionals qualified in 3D printing using DLP and LCD technologies.

1.4. Technical specifications and characteristics

1.4.1. Composition

Components
Methacrylate oligomers
Photoinitiator
Inorganic Filler
Pigments
Additives

2. STORAGE

Store in a dry place at a temperature between 35.6°F and 86°F. The product should not be exposed to sunlight.

3. INSTRUCTIONS FOR USE

3.1. Preparing the 3D Printing File

- Prepare the STL file in the software related to the 3D printer. Afterwards, it is suggested to initially print a calibration piece available on the Yllér website, using the indicated initial parameters.

3.2. Mixing the Resin

- Before adding the resin to the printer reservoir, place the bottle in a roller shaker for 30 minutes so that the inorganic particles of the filler are well distributed.

- After bringing the resin to room temperature in a water bath, keep it

within the temperature range of 82.4°F to 89.6°F to facilitate 3D printing.

- Add the appropriate amount of resin to the reservoir, depending on its current level. After printing is complete, return the resin to the bottle.

- **Note:** For the 50g version, add the total amount directly to the reservoir to make the device.

3.3. Post-processing

- Once the 3D printing process is complete, remove the printed devices from the platform. and carefully remove the supports. Then use one of the methods described below to remove excess resin from the device.

- Option 1: Place the device in a and add Cosmos Plus Composite Cleaning Solution until it covers the device and leave it for 1 minute, then brush it with a toothbrush for 1 minute. This process should be repeated 3 times, always changing the Cosmos Plus Composite Cleaning Solution. Afterwards, dry the device with compressed air.

- Option 2: Place the device in a beaker and add Cosmos Plus Composite Cleaning Solution until it covers the device and leave it for 3 minutes in the ultrasound, then remove it from the ultrasound and place it in a new Cosmos Plus Composite Cleaning Solution for 3 minutes at rest. Afterwards, dry the device with compressed air.

- **Anycubic post-curing (Wash & Cure 2.0 Model):** After drying, post-cure the device in the UV light chamber for 60 minutes. This is 30 minutes for each side.

- **Anycubic post-curing (Wash & Cure 3.0 Model):** After drying, post-cure the device in the UV light chamber for 60 minutes, with flexicure turned on. This is 30 minutes for each side.

- **Spectra Post-Curing (Model UV 01):** After drying, post-cure the

device in the UV light chamber with all LEDs on for 10 minutes at 73.4°F.

- **Pcure post-curing:** After drying, post-cure the device in the UV light chamber using the 'Upper' + 'Lower' 100% mode for 20 minutes.

- **UV Camera 72w or 36w post-curing:** After drying, post-cure the device in the UV light camera for 60 minutes. This is 30 minutes on each side.

- **Pionext post-curing (Model UV 01):** After drying, post-cure the device in the UV light camera, using all LEDs on for 10 minutes at 73.4°F.

- **Note:** Use the equipment indicated for post-curing or equipment with equivalent specifications.

- Finish: Polish to remove the support marks and then carry out the make-up process, according to the manufacturer's instructions.

- Finally, immerse the device in 70° alcohol for 10 minutes.

3.4. Cementation

- Carry out the protocol for preparing the prosthetic pieces.

- Apply Sylano (Yllor) according to the Instructions for Use.

- Prepare the teeth structure to receive the prosthetic piece and carry out the adhesive protocol. For best performance, we recommend using Astra Universal (Yllor). Follow the Instructions for Use.

- Proceed to cement the piece using the Ycem 4All Dual Resin Cement (Yllor) according to the Instructions for Use.

4. PRECAUTIONS, RESTRICTIONS AND WARNINGS

4.1. Precautions

- When handling this product, it is recommended to use Personal Protective Equipment, such as protective gloves, protective clothing, and eye protection.

- Make sure you carry out all the procedures mentioned in the instructions. Inefficient completion of any post-impression step may leave

traces of uncured liquid material that may alter the planning conditions of the device.

- Improper use of the printing equipment may compromise the physical and mechanical properties of the final device produced from this material. Follow the printer's Instructions for Use provided by the manufacturer and check the need for maintenance and calibration of the equipment.

- Keep the container closed when not in use.

- Keep this product out of the reach of children and animals.

- Do not use the product if its packaging is damaged.

- Do not use the product if it has expired.

- It is the responsibility of the professionals responsible to use the products and equipment in accordance with their Instructions for Use.

- Do not use the product without knowledge of the appropriate techniques. Always carry out the correct procedures in the right conditions and with the right equipment.

- Follow the Instructions for Use provided by the manufacturers of the software used for device planning.

- In the event of an accident involving this product, inform the manufacturer using the communication channels available on the website.

- Wash hands and other exposed areas thoroughly after handling.

- Avoid release into the environment.

4.2. Restrictions and warnings

The use of this product is contraindicated for people with signs of allergy or irritation to the chemical components of the material and/or methacrylates.

5. INDICATION

Product indicated for making veneers.

6. VALIDITY

Indicated on the packaging.

7. PRODUCT REGISTRATION

ANVISA: 81009359003.

8. DISPOSING OF MATERIALS

To dispose of materials, solid waste or remnants of the product, it is recommended to consult and comply with current legislation. Dispose of liquid contents and container at a special or hazardous waste collection point in accordance with local, regional, national and/or international regulations, at an authorized centre. If there is no legislation in force, pack them in a sharps collector and dispose of them in hospital waste.

Yller Biomaterials S/A is not responsible for damages, of any nature, caused by incorrect use or use not provided for in the Instructions for Use.

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